

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |   |
|--------------------------|---|
| <b>PRODUCT NAME</b>      | Recombinant Mouse CD268/TNFRSF13C Protein, N-His-SUMO |
| <b>PRODUCT REFERENCE</b> | ARO-P11064  |
| <b>ORIGIN SPECIES</b>    | Mus musculus (Mouse)                                  |
| <b>SOURCE</b>            |   |
| <b>FRAGMENT</b>          | Thr25-Ala71   |
| <b>TAG</b>               | N-His-SUMO Tag  |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>EXPRESSION SYSTEM</b>       | Prokaryotic expression   |
| <b>PRODUCTION HOST</b>         | <i>Escherichia coli (E. coli)</i>  |
| <b>MOLECULAR WEIGHT</b>        | 17.37 kDa  |
| <b>FORM</b>                    | Lyophilized  |
| <b>BUFFER</b>                  | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.   |
| <b>STABILITY &amp; STORAGE</b> | Store at 4°C for short term while lyophilized.<br>Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |