

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--------------------------------------------|
| PRODUCT NAME | Recombinant Mouse CD28 Protein, N-His-SUMO |
| PRODUCT REFERENCE | ARO-P11085 |
| ORIGIN SPECIES | Mus musculus (Mouse) |
| SOURCE | |
| FRAGMENT | Lys21-Lys137 |
| TAG | N-His-SUMO Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 25.91 kDa |
| FORM | Lyophilized |
| BUFFER | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |