



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Recombinant Human MBL2/MBP-C Protein, N-His-SUMO |
| PRODUCT REFERENCE | ARO-P11323 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| SOURCE | |
| FRAGMENT | Gly131-Ile248 |
| TAG | N-His-SUMO Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 25.30 kDa |
| FORM | Lyophilized |
| BUFFER | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY