

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant European house dust mite Der p 1/DERP1 Protein, N-His
PRODUCT REFERENCE	ARO-P11497
ORIGIN SPECIES	Dermatophagoides pteronyssinus (European house dust mite)
SOURCE	
FRAGMENT	Arg19-Leu320
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	36.57 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.