

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant JEV NS1/Non-structural protein 1 Protein, N-Trx-His & C-His
PRODUCT REFERENCE	ARO-P11533
ORIGIN SPECIES	Japanese encephalitis virus (strain Jaoars982) (JEV)
SOURCE	
FRAGMENT	Asp795-Ala1146
TAG	N-Trx-His & C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	58.42 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.