

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Recombinant Vibrio cholerae serotype O1 ctxB/Cholera Toxin Subunit B Protein, N-His-SUMO
<b>PRODUCT REFERENCE</b>	ARO-P11536
<b>ORIGIN SPECIES</b>	Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)
<b>SOURCE</b>	
<b>FRAGMENT</b>	Thr22-Asn124
<b>TAG</b>	N-His-SUMO Tag

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>PRODUCTION HOST</b>	<i>Escherichia coli</i> ( <i>E. coli</i> )
<b>MOLECULAR WEIGHT</b>	23.99 kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.