



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Human PICK1 Protein, N-His-SUMO & C-Strep
PRODUCT REFERENCE	ARO-P11839
ORIGIN SPECIES	Homo sapiens (Human)
SOURCE	
FRAGMENT	Gly15-Leu114
TAG	N-His-SUMO & C-Strep Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	23.98 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY