

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Recombinant Human HTR7 Protein, N-His-SUMO |
| PRODUCT REFERENCE | ARO-P11964 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| SOURCE | |
| FRAGMENT | Met258-Ala325 |
| TAG | N-His-SUMO Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 20.26 kDa |
| FORM | Lyophilized |
| BUFFER | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |