

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Recombinant Human NFIL3 Protein, N-His-SUMO & C-Strep |
| PRODUCT REFERENCE | ARO-P12277 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| SOURCE | |
| FRAGMENT | Glu65-Ser161 |
| TAG | N-His-SUMO & C-Strep Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 24.83 kDa |
| FORM | Lyophilized |
| BUFFER | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |