



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Recombinant EV71 VP1/Capsid protein VP1, N-His

PRODUCT REFERENCE ARO-P12706

ORIGIN SPECIES Human enterovirus 71 (EV71) (EV-71)

SOURCE

FRAGMENT Gly566-Leu862

TAG N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM Prokaryotic expression

PRODUCTION HOST *Escherichia coli* (*E. coli*)

MOLECULAR WEIGHT 35.20 kDa

FORM Lyophilized

BUFFER Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

STABILITY & STORAGE Store at 4°C for short term while lyophilized.
Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY