

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant DENV-2 Envelope protein E, N-His
PRODUCT REFERENCE	ARO-P12716
ORIGIN SPECIES	Dengue virus type 2 (strain Thailand/16681/1984) (DENV-2)
SOURCE	
FRAGMENT	Met281-Ala722
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	77.29 kDa
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.