



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Treponema pallidum/TP TPP15 + tmpA + TPP47 Protein, C-His
PRODUCT REFERENCE	ARO-P12738
ORIGIN SPECIES	Treponema pallidum (strain Nichols)
SOURCE	
FRAGMENT	TPP15[Ser19-Arg141] + tmpA[Gly20-Arg345] + TPP47[Gly21-Gln434]
TAG	C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	105 kDa
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY