

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	EndoS
PRODUCT REFERENCE	ARO-P12779
ORIGIN SPECIES	<i>Streptococcus pyogenes</i>
SOURCE	
FRAGMENT	EndoS is cloned from <i>Streptococcus pyogenes</i> and expressed in <i>E.coli</i> .
TAG	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli</i> (<i>E. coli</i>)
MOLECULAR WEIGHT	109.35 kDa
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.