

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant MDV Coat protein, N-His & N-SUMO
PRODUCT REFERENCE	ARO-P12847
ORIGIN SPECIES	MDV
SOURCE	
FRAGMENT	Val2-Leu172
TAG	N-His & N-SUMO Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli</i> (<i>E. coli</i>)
MOLECULAR WEIGHT	31.62 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.