

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Langya virus/LayV Polymerase Protein, N-His
PRODUCT REFERENCE	ARO-P15645
ORIGIN SPECIES	Langya virus (LayV)
SOURCE	UUV47207.1
FRAGMENT	Met1-Leu912
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	107.51 kDa
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.