



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	SA/Streptavidin Protein, HRP
PRODUCT REFERENCE	ARO-P15676
ORIGIN SPECIES	Streptomyces avidinii
SOURCE	P22629
FRAGMENT	Commercial secret.
TAG	HRP Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
PRODUCTION HOST	
MOLECULAR WEIGHT	
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY