



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Recombinant Human IgG3 Fc Fragment

PRODUCT REFERENCE ARO-P15944

ORIGIN SPECIES Homo sapiens (Human)

SOURCE P01860

FRAGMENT Glu99-Lys377

TAG

PRODUCT INFORMATION

EXPRESSION SYSTEM Eukaryotic expression

PRODUCTION HOST *Mammalian cells*

MOLECULAR WEIGHT 33.15 kDa

FORM Lyophilized

BUFFER Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

STABILITY & STORAGE Store at 4°C for short term while lyophilized.
Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY