

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant HAdV-2 E1B-19K Protein, N-GST & C-His
PRODUCT REFERENCE	ARO-P16501
ORIGIN SPECIES	Human adenovirus C serotype 2 (HAdV-2) (Human adenovirus 2)
SOURCE	P03247
FRAGMENT	Ser35-His108
TAG	N-GST & C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli</i> (<i>E. coli</i>)
MOLECULAR WEIGHT	36.67 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.