

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Recombinant KSI tag Control Protein (N-His)

PRODUCT REFERENCE ARO-P16749

ORIGIN SPECIES Pseudomonas testosteroni

SOURCE P00947

FRAGMENT Met1-Ala123

TAG

PRODUCT INFORMATION

EXPRESSION SYSTEM Prokaryotic expression

PRODUCTION HOST *Escherichia coli (E. coli)*

MOLECULAR WEIGHT 18.04 kDa

FORM Lyophilized

BUFFER Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

STABILITY & STORAGE Store at 4°C for short term while lyophilized.
Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Avoid freeze/thaw cycles, aliquot the product according to your needs.