



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant SARS-CoV-2 S - Trimer (RRAR-GSAS Protein, KV-PP), C-His
PRODUCT REFERENCE	ARO-P16896
ORIGIN SPECIES	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
SOURCE	P0DTC2
FRAGMENT	Gln14-Gln1208
TAG	C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
PRODUCTION HOST	<i>Mammalian cells</i>
MOLECULAR WEIGHT	139.82 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY