

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant SARS-CoV-2 ORF7a/Protein U122 Protein, N-His
PRODUCT REFERENCE	ARO-P16917
ORIGIN SPECIES	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
SOURCE	P0DTC7
FRAGMENT	Glu16-Tyr97
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	11.93 kDa
FORM	Lyophilized
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.