

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant HBV C/Capsid protein Protein, C-His
PRODUCT REFERENCE	ARO-P16931
ORIGIN SPECIES	Hepatitis B virus
SOURCE	QGD14736.1
FRAGMENT	Met1-Val149
TAG	C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli</i> (<i>E. coli</i>)
MOLECULAR WEIGHT	17.92 kDa
FORM	Lyophilized
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.