



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Endo S2/Endo-beta-N-acetylglucosaminidase
PRODUCT REFERENCE	ARO-P16934
ORIGIN SPECIES	Streptococcus pyogenes
SOURCE	
FRAGMENT	Endo S2/Endo-beta-N-acetylglucosaminidase is cloned from Streptococcus pyogenes and expressed in E.coli.
TAG	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	114.67 kDa
FORM	Lyophilized
BUFFER	Lyophilized from PBS pH7.4, 1mM EDTA, 4%trehalose, 1% mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY