

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Human SLC6A5 Protein, N-GST & C-His
PRODUCT REFERENCE	ARO-P17092
ORIGIN SPECIES	Homo sapiens (Human)
SOURCE	
FRAGMENT	Cys292-Arg393
TAG	N-GST & C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli</i> (<i>E. coli</i>)
MOLECULAR WEIGHT	39.81 kDa
FORM	Lyophilized
BUFFER	
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.