

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Recombinant Human rotavirus A Outer capsid protein VP4 Protein, N-His
PRODUCT REFERENCE	ARO-P17772
ORIGIN SPECIES	Human rotavirus A (RV-A)
SOURCE	
FRAGMENT	Ala248-Tyr480
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	28.64 kDa
FORM	Lyophilized
BUFFER	
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.