



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Mouse CASP9/Caspase 9 Recombinant Protein

PRODUCT REFERENCE ARO-P18431

ORIGIN SPECIES Mus musculus (Mouse)

SOURCE Q8C3Q9

FRAGMENT Asp196-Ser454

TAG

PRODUCT INFORMATION

EXPRESSION SYSTEM Prokaryotic expression

PRODUCTION HOST *Escherichia coli (E. coli)*

MOLECULAR WEIGHT

FORM Lyophilized

BUFFER PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

STABILITY & STORAGE Store at 4°C for short term while lyophilized.
Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY