

## PRODUCT DATASHEET

### OVERVIEW

---

|                   |                                                    |
|-------------------|----------------------------------------------------|
| PRODUCT NAME      | VZV/HHV-3 gH&gL Recombinant Protein, C-Strep&C-His |
| PRODUCT REFERENCE | ARO-P18739                                         |
| ORIGIN SPECIES    | HHV-3                                              |
| SOURCE            | Q775J3&Q9J3N1                                      |
| FRAGMENT          | Asn18-Arg795&Leu23-Gln160                          |
| TAG               | C-Strep&C-His Tag                                  |

### PRODUCT INFORMATION

---

|                     |                                                                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EXPRESSION SYSTEM   | Eukaryotic expression                                                                                                                                                                                                                        |
| PRODUCTION HOST     | <i>Mammalian Cells</i>                                                                                                                                                                                                                       |
| MOLECULAR WEIGHT    |                                                                                                                                                                                                                                              |
| FORM                | Lyophilized                                                                                                                                                                                                                                  |
| BUFFER              | PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                                                                                                            |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized.<br>Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |