

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Mouse CD11a/ITGAL Recombinant Protein, C-His
PRODUCT REFERENCE	ARO-P19191
ORIGIN SPECIES	Mus musculus (Mouse)
SOURCE	P24063
FRAGMENT	Lys151-Thr335
TAG	C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
PRODUCTION HOST	<i>Mammalian Cells</i>
MOLECULAR WEIGHT	24.06 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.