

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                       |
|--------------------------|---------------------------------------|
| <b>PRODUCT NAME</b>      | Human SIX3 Recombinant Protein, N-His |
| <b>PRODUCT REFERENCE</b> | ARO-P19765                            |
| <b>ORIGIN SPECIES</b>    | <i>Homo sapiens (Human)</i>           |
| <b>SOURCE</b>            | O95343                                |
| <b>FRAGMENT</b>          | Ser88-Gly276(Arg194Cys)               |
| <b>TAG</b>               | N-His Tag                             |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                                                                                                              |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>EXPRESSION SYSTEM</b>       | Prokaryotic expression                                                                                                                                                                                                                       |
| <b>PRODUCTION HOST</b>         | <i>Escherichia coli (E. coli)</i>                                                                                                                                                                                                            |
| <b>MOLECULAR WEIGHT</b>        | 24.26 kDa                                                                                                                                                                                                                                    |
| <b>FORM</b>                    | Lyophilized                                                                                                                                                                                                                                  |
| <b>BUFFER</b>                  | PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                                                                                                            |
| <b>STABILITY &amp; STORAGE</b> | Store at 4°C for short term while lyophilized.<br>Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |