



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human SUCNR1 Recombinant Protein, N-His-KSI
PRODUCT REFERENCE	ARO-P19922
ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
SOURCE	Q9BXA5
FRAGMENT	Leu2-Gly31&Arg81-His103&Pro159-Leu185&Tyr207-Pro230&Arg252-Arg281&Arg303-Lys334
TAG	N-His-KSI Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	35.81 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY