



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human BRS3 Recombinant Protein, N-His-KSI
<b>PRODUCT REFERENCE</b>	ARO-P19936
<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>SOURCE</b>	P32247
<b>FRAGMENT</b>	Ala2-Ser41&Asp104-Lys121&Phe185-Leu220&His294-Ile313&Ser334-Phe399
<b>TAG</b>	N-His-KSI Tag

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>PRODUCTION HOST</b>	<i>Escherichia coli (E. coli)</i>
<b>MOLECULAR WEIGHT</b>	36.26 kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**