



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|-----------------------------------------|
| PRODUCT NAME | Human ZNF124 Recombinant Protein, N-His |
| PRODUCT REFERENCE | ARO-P20045 |
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| SOURCE | Q15973 |
| FRAGMENT | Met1-Met351 |
| TAG | N-His Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 42.38 kDa |
| FORM | Lyophilized |
| BUFFER | PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY