

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human SLC5A6 Recombinant Protein, N-His-KSI
PRODUCT REFERENCE	ARO-P20178
ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
SOURCE	Q9Y289
FRAGMENT	Tyr122-Thr140&Asn274-Ser293&Ala317-Pro352&Gly550-Leu636
TAG	N-His-KSI Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	34.35 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.