

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Human GPR37L1 Recombinant Protein, N-His-KSI |
| PRODUCT REFERENCE | ARO-P20214 |
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| SOURCE | O60883 |
| FRAGMENT | Ala26-Met134 |
| TAG | N-His-KSI Tag |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| PRODUCTION HOST | <i>Escherichia coli (E. coli)</i> |
| MOLECULAR WEIGHT | 27.68 kDa |
| FORM | Lyophilized |
| BUFFER | PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol. |
| STABILITY & STORAGE | Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs. |