

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human STAB1/Stabilin-1 Recombinant Protein, N-His
PRODUCT REFERENCE	ARO-P20529
ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
SOURCE	Q9NY15
FRAGMENT	Gly357-Gly492
TAG	N-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	17.68 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.