



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

products@proteogenix.fr

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT NAME</b>      | Human AKAP7 Recombinant Protein, N-GST & C-His |
| <b>PRODUCT REFERENCE</b> | ARO-P20683                                     |
| <b>ORIGIN SPECIES</b>    | <i>Homo sapiens (Human)</i>                    |
| <b>SOURCE</b>            | O43687   |
| <b>FRAGMENT</b>          | Ile16-Lys104                                   |
| <b>TAG</b>               | N-GST & C-His Tag                              |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>EXPRESSION SYSTEM</b>       | Prokaryotic expression   |
| <b>PRODUCTION HOST</b>         | <i>Escherichia coli (E. coli)</i>  |
| <b>MOLECULAR WEIGHT</b>        | 37.95 kDa  |
| <b>FORM</b>                    | Lyophilized  |
| <b>BUFFER</b>                  | PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.  |
| <b>STABILITY &amp; STORAGE</b> | Store at 4°C for short term while lyophilized.<br>Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**