



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD365/KIM-1/TIM-1/HAVCR1 Recombinant Protein, C-Fc
PRODUCT REFERENCE	ARO-P20785
ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
SOURCE	Q96D42
FRAGMENT	Ser21-Thr293
TAG	C-Fc Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
PRODUCTION HOST	<i>Mammalian Cells</i>
MOLECULAR WEIGHT	57.09 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY