

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Mouse CCL7/MCP-3 Recombinant Protein, N-GST & C-His
PRODUCT REFERENCE	ARO-P21074
ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
SOURCE	Q03366
FRAGMENT	Lys40-Pro97
TAG	N-GST & C-His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
PRODUCTION HOST	<i>Escherichia coli (E. coli)</i>
MOLECULAR WEIGHT	34.97 kDa
FORM	Lyophilized
BUFFER	PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	Store at 4°C for short term while lyophilized. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Avoid freeze/thaw cycles, aliquot the product according to your needs.