

Datasheet
Aflibercept Biosimilar - Anti-VEGF mAb - Research Grade
PX-TA1016

➤ **Overview**

Antibody name: Aflibercept Biosimilar - Anti-VEGF mAb - Research Grade
Isotype: hG1-Fc
Host species: Human
Production Host: Mammalian cells
Source: DrugBank DB08885
Aliases/Synonyms: Aflibercept

➤ **Properties**

Clonality: Monoclonal Antibody
Final buffer: PBS pH7.5
Purity: ≥90%
Molecular Weight: 115kDa

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.