

Datasheet

Tocilizumab Biosimilar - Anti-IL-6R mAb - Research Grade PX-TA1030

➤ Overview

| | |
|--------------------------|--|
| Antibody name: | Tocilizumab Biosimilar - Anti-IL-6R mAb - Research Grade |
| Isotype: | IgG1Kappa |
| Host species: | Humanized |
| Production Host: | Mammalian cells |
| Source: | IMGT8394 |
| Aliases/Synonyms: | Tocilizumab |

➤ Properties

| | |
|--------------------------|---------------------|
| Clonality: | Monoclonal Antibody |
| Final buffer: | PBS pH7.5 |
| Purity: | ≥90% |
| Molecular Weight: | 150kDa |

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.