

Datasheet

4F2Mab Biosimilar - Anti-ganglioside GD2/ SFRP1/SLC3A2 mAb - Research Grade PX-TA1073

➤ Overview

Antibody name: 4F2Mab Biosimilar - Anti-ganglioside GD2/ SFRP1/SLC3A2 mAb - Research Grade

Isotype: IgG2a-nd; IgG1-nd

Host species: Mus musculus, Humanized

Production Host: Mammalian cells

Aliases/Synonyms: 14G2a, 4F2

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity: >95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.