

Datasheet

Firivumab Biosimilar - Anti-hemagglutinin HA mAb - Research Grade PX-TA1361

➤ Overview

| | |
|--------------------------|---|
| Antibody name: | Firivumab Biosimilar - Anti-hemagglutinin HA mAb - Research Grade |
| Isotype: | IgG1-kappa |
| Host species: | Homo sapiens |
| Production Host: | Mammalian cells |
| Source: | CAS 1443004-15-6 |
| Aliases/Synonyms: | CT-P22 |

➤ Properties

| | |
|----------------------|---------------------|
| Clonality: | Monoclonal Antibody |
| Final buffer: | PBS buffer PH7.5 |
| Purity: | >95% |

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.