

## Datasheet

### Lodapolimab Biosimilar - Anti-PDCD1/PD-1/CD279 mAb - Research Grade PX-TA1593

#### ➤ Overview

**Antibody name:** Lodapolimab Biosimilar - Anti-PDCD1/PD-1/CD279 mAb - Research Grade

**Isotype:** IgG1-lambda

**Host species:** Homo sapiens

**Production Host:** Mammalian cells

**Source:** CAS 2118349-31-6

**Aliases/Synonyms:** LY3300054

#### ➤ Properties

**Clonality:** Monoclonal Antibody

**Final buffer:** PBS buffer PH7.5

**Purity:** >95%

## Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

## Note:

For research use only. Not suitable for clinical or therapeutic use.