

Datasheet

Demupitamab Biosimilar - Anti-EGFR mAb - Research Grade

PX-TA1655

➤ Overview

Antibody name: Demupitamab Biosimilar - Anti-EGFR mAb - Research Grade

Isotype: IgG1, Kappa

Host species: Homo sapiens

Production Host: XtenCHO

Source: CAS 2249926-74-5

Aliases/Synonyms: IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS EGFR (EPIDERMAL GROWTH FACTOR RECEPTOR, RECEPTOR TYROSINE-PROTEIN KINASE ERBB-1, ERBB1, HER1, HER-1, ERBB)), HOMO SAPIENS MONOCLONAL ANTIBODYIMMUNOGLOBULIN IGG1, ANTI-(HUMAN EPIDERMAL GROWTH FACTOR RECEPTOR) (HU

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity: >95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.