

**Datasheet**  
**Ezabentimab Biosimilar - Anti-PDCD1 mAb - Research Grade**  
**PX-TA1666**

➤ **Overview**

|                          |   |
|--------------------------|---|
| <b>Antibody name:</b>    | Ezabentimab Biosimilar - Anti-PDCD1 mAb - Research Grade                        |
| <b>Isotype:</b>          | IgG4,Kappa  |
| <b>Host species:</b>     | Humanized   |
| <b>Production Host:</b>  | XtenCHO   |
| <b>Source:</b>           | CAS 2249882-54-8  |
| <b>Aliases/Synonyms:</b> | ANTI-PD-1 MONOCLONAL ANTIBODY BI 754091,BI-754091,EZABENTIMAB,EZABENTIMAB [INN] |

➤ **Properties**

|                      |                     |
|----------------------|---------------------|
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>Final buffer:</b> | PBS buffer PH7.5    |
| <b>Purity:</b>       | >95%                |

**Storage advice:**

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

**Note:**

For research use only. Not suitable for clinical or therapeutic use.