

Datasheet

Lemzoparlimab Biosimilar - Anti-CD47 mAb - Research Grade PX-TA1677

➤ Overview

Antibody name:	Lemzoparlimab Biosimilar - Anti-CD47 mAb - Research Grade
Isotype:	IgG4,Kappa
Host species:	Humanized
Production Host:	XtenCHO
Source:	CAS 2377483-71-9
Aliases/Synonyms:	IMMUNOGLOBULIN G4(228-PROLINE), ANTI-(HUMAN CD47 ANTIGEN) (HUMAN MONOCLONAL TJ011133 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN MONOCLONAL TJ011133 .KAPPA.-CHAIN, DIMER,TJ011133, TJ-011133

➤ Properties

Clonality:	Monoclonal Antibody
Final buffer:	PBS buffer PH7.5
Purity:	>95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.