

Datasheet

Lorukafusp Biosimilar - Anti-Ganglioside GD2 mAb - Research Grade PX-TA1681

➤ Overview

Antibody name: Lorukafusp Biosimilar - Anti-Ganglioside GD2 mAb - Research Grade

Isotype: IgG1, Kappa

Host species: Chimeric

Production Host: XtenCHO

Source: CAS 2131168-99-3

Aliases/Synonyms: ANTI-GANGLIOSIDE-GD2-ANTIBODY-INTERLEUKIN-2-FUSION-PROTEIN, APN301, EMD 273063, HU14.18-IL2, IMMUNOGLOBULIN G1-KAPPA, ANTI-(HUMAN GANGLIOSIDE GD2) CHIMERIC MONOCLONAL ANTIBODY, FUSED WITH INTERLEUKIN 2, GLYCOFORM ALFA, APN-301, EMD-273063, HU 14.18IL2

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity: >95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.