

Datasheet

Ociperlimab Biosimilar - Anti-TIGIT mAb - Research Grade

PX-TA1693

➤ Overview

Antibody name: Ociperlimab Biosimilar - Anti-TIGIT mAb - Research Grade

Isotype: IgG1,Kappa

Host species: Humanized

Production Host: XtenCHO

Source: CAS 2342597-93-5

Aliases/Synonyms: BGB-A1217, IMMUNOGLOBULIN G1, ANTI-(HUMAN RECEPTOR TIGIT) (HUMAN-MUS MUSCULUS MONOCLONAL BGB-A1217 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL BGB-A1217 .KAPPA.-CHAIN, DIMER, IMMUNOGLOBULIN G1-KAPPA, ANTI-(HOMO SAPIENS TIGIT (T-CELL IMMU

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity: >95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.