

Datasheet

Sabatolimab Biosimilar - Anti-HAVCR2 mAb - Research Grade PX-TA1714

➤ Overview

Antibody name:	Sabatolimab Biosimilar - Anti-HAVCR2 mAb - Research Grade
Isotype:	IgG4,Kappa
Host species:	Humanized
Production Host:	XtenCHO
Source:	CAS 2252262-24-9
Aliases/Synonyms:	ANTI-TIM3 CHECKPOINT INHIBITOR MBG453, IMMUNOGLOBULIN G4 (226-PROLINE,DE-C-TERMINAL- LYSINE), ANTI-(HUMAN HEPATITIS VIRUS A CELLULAR RECEPTOR 2) (HUMAN-MUS MUSCULUS MONOCLONAL CLONE NVS260714 .GAMMA.4-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL CL

➤ Properties

Clonality:	Monoclonal Antibody
Final buffer:	PBS buffer PH7.5
Purity:	>95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.