

Datasheet

Upifitamab Biosimilar - Anti-SLC34A2 mAb - Research Grade

PX-TA1732

➤ Overview

Antibody name: Upifitamab Biosimilar - Anti-SLC34A2 mAb - Research Grade

Isotype: IgG1,Kappa

Host species: Humanized

Production Host: XtenCHO

Source: CAS 2254118-43-7

Aliases/Synonyms: IMMUNOGLOBULIN G1 (DE-450-LYSINE), ANTI-(HUMAN SODIUM-DEPENDENT PHOSPHATE TRANSPORT PROTEIN 2B) (HUMAN-MUS MUSCULUS MONOCLONAL XMT-1535 .GAMMA.1-CHAIN), DISULFIDE WITH HUMAN-MUS MUSCULUS MONOCLONAL XMT-1535 .KAPPA.-CHAIN, DIMER,

➤ Properties

Clonality: Monoclonal Antibody

Final buffer: PBS buffer PH7.5

Purity: >95%

Storage advice:

ProteoGenix did not perform stability or storage tests for your specific end products.

Here are some **standard storage advice**:

- Store at 4°C for short term; Freeze and store at -20°C or -80°C for long term
- Freezing should be first tested on small aliquot(s); Glycerol up to 50% may be added for cryoprotection. Purified protein/antibodies sometimes require additives for optimal freezing and storage, and/or a specific freezing process (e.g., -20°C, -80°C, or flash freezing in liquid nitrogen).
- Purified proteins/antibodies should ideally be stored as 0.5 to 2mg/ml stock solutions.
- Avoid freeze/thaw cycles: aliquot the products according to your needs.

Note:

For research use only. Not suitable for clinical or therapeutic use.